



Peanut Massage Ball

Designed to hit those hard-to-reach spots of the body, utilize the contoured Peanut Massage Ball to target areas including the neck, forearms, Achilles tendon, shoulders, calves and more. Suitable for reflexology and both home and therapeutic use, the eco-conscious Peanut Massage Ball's gentle spikes stimulates blood flow and promotes circulation while loosening muscle knots.

Includes 5 Exercises:

- Neck Massage
- Forearm Extensor Release
- Forearm Flexor Release
- Achilles Tendon Release
- Achilles Tendon Deep Release



merrithew™ | leaders en
mouvement conscient...

Neck Massage



Starting position: Lie facing up on a mat. Place Peanut Massage Ball just above the nape of the neck.



Nod head forward dropping the chin toward the chest.



Lift chin up toward the ceiling, extending the neck.

Forearm Extensor Release



Starting position: Kneel on a mat, hands down for support. Place Peanut Massage Ball under back of one forearm with palm up.



Roll the Massage Ball from the wrist to the elbow applying slight pressure.

Forearm Flexor Release



Starting position: Kneel on a mat, hands down for support. Place Peanut Massage Ball under front of one forearm with palm down.



Roll the Massage Ball from the wrist to the elbow applying slight pressure.

Achilles Tendon Release



Starting position: Sit on a mat in a comfortable position, hands behind torso supporting weight.

1. Bend knees, place one heel on Peanut Massage Ball.
2. Lengthen leg, applying pressure on Massage Ball and rolling it up the tendon.
3. Bend knee to return.
4. Increase the degree of stretch by keeping the ankle in dorsiflexion



WARRANTY, CARE AND USAGE FOR PEANUT MASSAGE BALL

Warranty

This product is sold without any warranties or guarantees of any kind. The manufacturer and distributor disclaim any liability, loss or damage caused by its use.

Care and Cleaning

- Wipe clean using a diluted mixture of water and mild soap
- Rinse well
- Ensure product is dry before using
- Never use corrosive cleaning products

Safety and Usage

Improper use of exercise equipment may cause serious bodily injury. To reduce risk, please read the following information carefully.

- Before starting any exercise program, consult a physician.
- Stop exercising if you experience chest pain, feel faint, have difficulty breathing, or experience muscular or skeletal discomfort.
- Perform exercises in a slow and controlled manner; avoid excessive tension in neck and shoulders.
- Do not use equipment if it appears worn, broken or damaged, and do not attempt to repair equipment yourself.
- Use equipment only as recommended.

Achilles Tendon Deep Release



Starting position: sit on a mat in a comfortable position, hands behind torso supporting weight.

1. Bend knees with one heel on Peanut Massage Ball, place the heel of the opposite foot on the toes.
2. Lengthen leg, applying pressure on Massage Ball and rolling it up the tendon.
3. Bend knee to return.
4. Increase the degree of stretch by keeping the ankle in dorsiflexion



2200 Yonge Street, Suite 500
Toronto, Ontario, Canada M4S 2C6
416.482.4050 | 1.800.910.0001

