Total Body Massage Stick




Muscle groups and body parts that can be targeted:

Hip Flexor
Quadriceps
Anterior Chest
Posterior Shoulder
Upper Arm
Latissimus Dorsi
Deltoids
Trapezius
Inside Forearm
Outside Forearm
Calf
Anterior Shin

Treating trigger points and achieving myofascial release are simple and accessible with the Total Body Massage Stick. Slow, steady pressure applied to 'knots' in a muscle allow it to relax thereby increasing the flow of blood to the site. The Total Body Massage Stick features a row of independently spinning cylinders that can localize the pressure in any part of the body for a more effective result.

Muscle groups and body parts that can be targeted:

To use the Massage Stick, grasp the stick with both hands and roll it firmly over the target muscle up to 20 times using short strokes. Avoid rolling over bony structures, staying on the muscle belly. The Massage Stick can be used directly on the skin or through a thin layer of clothing. Keep muscles as relaxed as possible during the application and apply only as much pressure as can be tolerated comfortably. **You should not feel pain.**

CARE, USAGE, AND WARRANTY FOR MASSAGE STICK

Care and Cleaning

If necessary, clean the Massage Stick with a combination of tea tree oil* and water. A mixture of mild soap and water can be used to remove more persistent dirt. Ensure cleaner used does not make surfaces slippery. Do not use corrosive cleaning products on the Massage Stick.

*Tea tree oil is a natural disinfectant available at most pharmacies or health food stores; mix according to directions on package.

Safety and Usage

Improper use of exercise equipment may cause serious bodily injury. To reduce risk, please read the following information carefully.

- The information provided in this exercise guide is for educational purposes only.
- The Massage Stick should only be used under the advice of a medical practitioner as tissue massage may not be recommended.
- Gentle pressure should be used to avoid damaging the muscle tissue.
- Periods of extended use may result in muscle soreness.
- Mild discomfort should be felt when rolling over trigger points – never pain.
- Do not use on skin that is irritated damaged or broken.
- Muscles should remain relaxed during use.
- Do not use this product if it appears worn, broken or damaged, and do not attempt to repair equipment yourself.

Warranty

This product is sold without any warranties or guarantees of any kind. The manufacturer and distributor disclaim any liability, loss or damage caused by its use.

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